

TauroLock™

Catalogue # TP-01

A. Description and Specifications

TauroLock™ contains anticoagulant and antimicrobial substances. It is to be used with devices for venous access e.g. a port or a catheter-based vascular access device. It is to be instilled in the device lumens between treatments in order to make the internal flow passages resistant to clot formation and hostile to bacterial and fungal growth. The solution must be withdrawn prior to initiating the next treatment. Active ingredients in TauroLock™ are (cyclo)-taurolidine and citrate (4%). Other components include water for injection and PVP. The pH is adjusted with citrate and/or sodium hydroxide. The product is sterile filter processed and supplied as a clear, sterile, non-pyrogenic solution.

Note:

For complete details of catheter-based venous vascular access products, consult the manufacturer's instructions or clinician's manual.

B. Indications

TauroLock™ is indicated for those patients who use a port or a silicone or polyurethane catheter-based device for their venous vascular access. TauroLock™ is intended to be used as a catheter lock solution by healthcare professionals. It is to be instilled into the device at the termination of a treatment and withdrawn prior to initiating subsequent treatments (see F4).

C. Contraindications

TauroLock™ is contraindicated for patients with a known allergy to citrate or (cyclo)-taurolidine, or when a patient is currently taking medication with known adverse interaction to citrate or (cyclo)-taurolidine.

D. Cautions

- As a consumable TauroLock™ is for single use only. Reuse creates a potential contamination risk for the patient.
- TauroLock™ is not for systemic injection. TauroLock™ must be used as a catheter lock solution as described in the access device's instruction for use. Failure to adhere to these instructions may result in inadvertent systemic injection of the solution. Once instilled into the catheter the solution must not be used again after aspiration.
- The vial* is a multi dose container. Once punctured it must be used within 48 hours. The ampoule is for single dose only due to potential risk of contamination. Multi use is contraindicated due to potential risk of contamination.
- In the event that access device patency is compromised, follow institutional protocol for restoring flow. Note: In case of patency problems other TauroLock™ variants such as TauroLock™-HEP100, TauroLock™-HEP500 or TauroLock™-U25.000*, which contain additional anticoagulant or fibrinolytic agents, are available.
- The exact filling volume of the access device has to be strictly respected with infants and children less than two years old, due to the active ingredient of citrate.
- In access devices which have previously been locked with non-antimicrobial lock solutions (e.g. with heparin, low concentrated citrate or saline) prior to the use of TauroLock™, viable organisms and endotoxins may be released from the biofilm. In this situation the lock solution must be aspirated prior to the next treatment to prevent very rare anaphylactic reactions which are not attributable to the active ingredients.
- The concentration of the antimicrobial compound is near to saturation. If not stored or transported according to the instructions under section H, precipitation can occur in the product. Do not use such a precipitated product.
- Blood drawn from catheters locked with TauroLock™ shall not be used for measurement of blood parameters (due to potential falsification).

E. Adverse Effects

Assessment of adverse effects is based on the following definitions of incidence:

| Very common | Common | Uncommon | Rare | Very rare | Not known |
|-------------|------------------|---------------------|-------------------------|-------------|---|
| ≥ 1/10 | ≥ 1/100 - < 1/10 | ≥ 1/1.000 - < 1/100 | ≥ 1/100.000 - < 1/1.000 | < 1/100.000 | cannot be estimated from the available data |

The following undesired effects may occur: Anaphylaxis (very rare); mild Hypocalcemia (common). There are no known risks associated with concomitant systemic antibiotic therapy or exposure to magnetic fields.

F. Instillation of TauroLock™

Follow the manufacturer's instructions that accompany the particular venous vascular access product utilized. Specific catheter lock volumes are associated with each device.

- Flush the device with 10 mL of saline.
- Withdraw TauroLock™ from the container using an appropriate syringe.
- Instill TauroLock™ slowly (not more than 1 mL per second, infants and children less than two years of age not more than 1 mL per 5 seconds) into the access device in a quantity sufficient to fill the lumen completely. **Consult the manufacturer's instructions for the specific fill volume or specify fill volume during implantation. The volume has to be strictly respected.** TauroLock™ will remain inside the access device until the next treatment (up to a maximum of 30 days).
- Prior to the next treatment, TauroLock™ must be aspirated and discarded in accordance with the institution's policy for infectious waste disposal. In case TauroLock™ cannot be aspirated slow flushing (not more than 1 mL per 3 seconds) is possible. If in exceptional cases the healthcare professional decides that aspiration of TauroLock™ is not appropriate (blood in the catheter, e. g. in parenteral nutrition) slow flushing of TauroLock™ is clinically tolerable. Do not flush TauroLock™ in infants and children less than two years of age, even in exceptional cases.
- Flush the device with 10 mL of saline.

G. Pregnancy and Breastfeeding

No data are available for pregnant and breastfeeding women. For safety reasons TauroLock™ should not be used during pregnancy and breastfeeding.

H. Storage and shipment

TauroLock™ must be stored at a temperature of 15 to 30°C and must not be shipped at freezing temperature. Do not freeze.

I. Packaging configuration

The following packaging configurations are available for TauroLock™: 10 x 3 mL TauroLock™ ampoules (single dose container). 10 x 5 mL TauroLock™ ampoules (single dose container). 100 x 10 mL TauroLock™ vials* (multi dose container).

*Not available in Australia.

State: 17.05.2019

 TauroPharm GmbH · August-Bebel-Straße 51 · D-97297 Waldbüttelbrunn · Germany
Tel: +49 931 304 299 0 · Fax: +49 931 304 299 29



Sterile, aseptic fill.



Read instruction for use.



Single use. The ampoule is a single dose and the vial* a multi dose container.



Do not use when package is damaged.



0123

CE acc. MDD 93/42/EEC, notified body:
TÜV SÜD PRODUCT SERVICE GmbH.